

# SOP FOR ELECTROTRANSFER OF SDS-POLYACRYLAMIDE GELS

## **TRANSFER BUFFER**

CAPS buffer is preferred for N-terminal sequence analysis because the background is significantly lower and there is no glycine present to contaminate the first cycle of automated Edman degradation.

**CAPS: 10 mM CAPS, 10% methanol, pH 11.0. Degas before use.**  
(CAPS = 3-[cyclohexylamino]-1-propane-sulfonic acid); Sigma C-6070)

**To prepare 1 liter of a 10X concentrated stock solution:**

- ☐ Place 900 ml Milli-Q water in a 1 liter beaker; insert stir bar
- ☐ Dissolve 22.13 gm CAPS (chemical shelf, RT) in water
- ☐ Titrate dropwise with 10 M NaOH to final pH
- ☐ Transfer to a 1000 ml graduated cylinder
- ☐ Add Milli-Q water to 1000 ml
- ☐ Store at 4C

**To prepare 1 liter of working solution:**

- ☐ Measure 800 ml of Milli-Q water into a 1000 ml graduated cylinder
- ☐ Add 100 ml of 10X stock CAPS buffer
- ☐ Add 100 ml 100% HPLC grade methanol
- ☐ Mix well
- ☐ Store at RT

## **GEL AND MEMBRANE EQUILIBRATION**

☐ Cut PVDF membrane to the same size as the gel. Pre-wet the PVDF by immersion in 100% methanol for 30 sec, followed by equilibration in transfer buffer for 5-15 min.

☐ Equilibrate the gel in transfer for 15-30 min (depending upon gel thickness).

*This is actually an important step. Many scientists do not equilibrate their gel. This leads to a drastic reduction in blotting efficiency.*

## **BLOTTING**

**IT IS CRITICALLY IMPORTANT TO NOT CAPTURE AIR BUBBLES BETWEEN THE MEMBRANE AND THE GEL. THIS WILL BLOCK TRANSFER OF PROTEINS TO THE MEMBRANE.**

### **MATERIALS:**

- ☐ One Large staining tray
- ☐ One Small staining tray with CAPS buffer
- ☐ Equilibrated PVDF membrane
- ☐ Equilibrated gel
- ☐ 10 x 10 cm gel blot paper-2 pieces
- ☐ Mini-gel blotting sandwich cassette materials (frame, 2 sponges)
- ☐ 1X CAPS buffer

### **PROCEDURE:**

- ☐ Pour a small amount of CAPS buffer into a small staining tray
- ☐ Place one sponge in tray
- ☐ Place two sheets of gel blot paper in tray
- ☐ Place another sponge in tray
- ☐ Make certain that the transfer tank has CAPS buffer
- ☐ Place the blotting sandwich frame into the large staining gray with the **black side down**
- ☐ Place one sponge on the frame
- ☐ Pour enough CAPS buffer into the tray to just cover the frame and sponge
- ☐ Place a piece of soaked gel blot paper on the sponge
- ☐ Place gel on paper
- ☐ Place equilibrated PVDF membrane on the gel
- ☐ Place another piece of soaked gel blot paper on the membrane
- ☐ Place a soaked sponge on the paper
- ☐ Close and snap shut the sandwich cassette
- ☐ Transfer immediately to the transfer tank

Transfer settings: **ALL GELS TRANSFERRED AT 100 mamps (constant)**

10% acrylamide, 0.75 mm spacer	1 hour
15% acrylamide, 0.75 mm spacer	1.5 hour
10-15% acrylamide, 1 mm spacer	2 hours
10-15% acrylamide, 1.5 mm spacer	3 hours