

# SOP FOR NON-FIXING SILVER STAIN FOR SDS POLYACRYLAMIDE GELS (BLUM METHOD)

(Ref: Blum, H., Beier, H. and Gross, H.J. (1986). Improved silver staining of plant proteins, RNA and DNA in polyacrylamide gels. Electrophoresis 8: 93-99.)

## SOLUTIONS

**Fixative:** Acetic acid 10% (v/v), ethanol 30% (v/v)

*Measure 100 ml Glacial Acetic Acid  
Measure 300 ml ethanol  
To 1000 ml with MilliQ water*

**Rinse:** 20% (v/v) ethanol

**Sensitizer:** 0.02% (w/v) sodium thiosulfate

*Weigh 0.2 gm sodium thiosulfate  
Dissolve in 1000 ml MilliQ water*

**Silver nitrate:** 0.2% (w/v) silver nitrate

*Measure 2.0 gm silver nitrate  
Dissolve in 1000 ml MilliQ water*

**Developer:** MUST BE PREPARED FRESH

3% (w/v) sodium carbonate  
0.025% (v/v) formaldehyde  
sodium thiosulfate 10 mgm/L

*Measure 15 gm sodium carbonate  
Measure 25 ml of stock sodium thiosulfate  
Measure 125 ul 37% formaldehyde  
To 500 ml with MilliQ water*

**Stop:** Tris, acetic acid

*Measure 50 gm Tris base  
Measure 25 ml Glacial acetic acid  
To 1000 ml with MilliQ water*

## PROCEDURE

1. Soak the gel in Fixative for at least one hour; change solution for a minimum of another hour; overnight is ok.

2. Rinse gel in Rinse Solution for 20 min.
3. Rinse gel in MilliQ water for 10 min (minigel) or 20 min (13 x 16cm gel).
4. Soak gel in Sensitizer Solution for 1 min.
5. Rinse gel in MilliQ water for 3 x 20 seconds.
6. Soak gel in Silver Nitrate Solution for 45 min. **(MAKE DEVELOPER NOW)**
7. Rinse gel with MilliQ water 5-10 seconds. **THIS IS A CRITICAL TIME. IT MEANS SECONDS, NOT MINUTES.**
8. Soak gel in Developer Solution until bands are adequate (typically 8-10 minutes)
9. Soak gel in Stop Solution for minimum 15 minutes, then keep in water.

**THE GEL STAINING PROCESS IS COMPLETE. IT DOES NOT REQUIRE DESTAINING LIKE A COOMASSIE STAINED GEL. DESTAIN ONLY SLICES THAT WILL BE DIGESTED.**

#### **DESTAINING SILVER STAINED GEL SLICES (for Blum silver stained gels only)**

Solutions needed: 30 mM potassium ferricyanide; 98 mg/10 ml water  
(make fresh) 100 mM sodium thiosulfate; 248 mg/10 ml water

#### Procedure:

mix solutions together 1:1 (this is working solution)  
Pipet 50  $\mu$ l of working solution onto gel slice(s)  
Incubate until slice is destained (several minutes); stop with water; discard  
Cover with 250  $\mu$ l of 100mM ammonium bicarb, 20 min; discard